

SEEDBURO PRECISION DIVIDER

Seedburo offers the traditional positive action divider that handles seed, grain and feed. The Precision Divider has repeatedly demonstrated its ability to precisely divide larger samples into working size samples in a matter of seconds. Ideal for smaller seeded crops.

The dividing and mixing are accomplished by the centrifugal action of a motor driven revolving neoprene disc positioned under the hopper. A smooth inside chrome plated finish virtually eliminates any of the sample sticking to the inner surfaces. Both the hopper and cylinder top lift off to provide easy access to the interior of the dividing chamber. Actual dims. 13" L x 13" W x 19-1/4" H, Net wt. 41 lbs., Ship wt. 46 lbs., Ship dims. 15" L x 15" W x 22" H.

Precision divider, 1/12HP, 115V, 60Hz, 1PH motor No. 106
Transformer for 220V operation No. K22



LABORATORY ROTARY DIVIDERS

The Gamet Rotary Dividers are fast and accurate to within 1/10th of 1% of the total sample by weight. Minimum sample sizes are 1000 grams. These dividers can significantly reduce the length of time it takes one employee to handle and divide larger samples. These dividers eliminate the need to divide samples many times as with the traditional gravity type dividers.

Our models 4203 and 4304 have 10 rotary divisions. When placing an order for the 4203 unit, **please specify the divisions you require.** Any combination of three numbers adding up to 10, for example a division of 3-4-3 or 5-3-2. The model 4204 again has 10 rotary divisions, but divides the sample into 4 portions. Any combination of 4 numbers adding up to 10 will tell us the divisions you require. The model 4503 has 12 divisions and divides into 3 portions. Therefore, any combinations of 3 numbers that add up to 12 will tell us the divisions you require.

Model 4203 dims. 22" L x 22½" W x 33" H, Net wt. 75 lbs. Model 4204 dims. 22" L x 22½" W x 33" H. Net wt. 85 lbs. Export packed: Ship wt. 150 lbs., Ship dims. 26" L x 26" W x 39" H

3-Way Laboratory Rotary Divider, 1/4HP, 115V, 50/60Hz, 1PH No. R4203
3-Way Laboratory Rotary Divider, 1/4HP, 230V, 50/60Hz, 1PH No. R4203/A
4-Way Laboratory Rotary Divider, 1/4HP, 115V, 50/60Hz, 1PH No. R4204
4-Way Laboratory Rotary Divider, 1/4HP, 230V, 50/60Hz, 1PH No. R4204/A
Export Packaging No. R4200/EX



RIFFLE DIVIDER

The Riffle Divider is used for halving dry materials. As the product is poured into the hopper, it is divided into two equal portions by a series of chutes which discharge alternately in opposite directions. Prices include stainless steel dividers, cold rolled steel pans, scoop and cleaning brush. The R75 is enclosed and all other models are open style as shown in the photograph. Extra scoops and brushes are available from Seedburo, please contact us for prices.

Number of Chutes	Width of Chutes	Pans Furnished	Hopper Size	Extra Pans	Model Number
14	3/8"	4	8" x 6-3/4"	R62P	R62
14	1 / 2"	4	11" x 6-3/4"	R64P	R64
14	3 / 4"	4	14-3/4" x 6-3/4"	R66P	R66
24	3 / 8"	2	4-1/2" x 11"	R75P	R75
12	1 / 2 "	4	11" x 8-1/2"	R80P	R80
12	3 / 4"	4	14" x 8-1/2"	R85P	R85
16	1"	4	20" x 9"	R87P	R87
8	2"	4	19-1/2" x 9-1/2"	R89P	R89
10	1-1/2"	4	20-1/2" x 9"	R90P	R90
8	2-1/2"	4	24" x 9"	R92P	R92*

*Pans with this unit are stainless steel.

