



SEEDBURO COUNTING AND PLACEMENT TRAYS

The Seedburo Counting and Placement Trays are constructed of durable acrylic material with brass screws. Each 2-piece tray features a spring-loaded encasement that allows for quick, uniform discharge of seed. Ideal for placement on germination pads or for quick seed counts. These trays are washdown safe to avoid cross contamination or treatment build-up. Standard trays are pre-drilled with 100 openings and are available in three frame sizes. **Custom frame size, hole size and hole count available, contact us for quote.**



Model	Frame Size	Commodity	Commonly Used w/Tray Size
No. 99B	9" x 9"	Soybean	19.5" x 18.75"
No. 99C	9" x 9"	Corn	19.5" x 18.75"
No. 912B	9" x 12"	Soybean	16" x 24" Food Tray
No. 912C	9" x 12"	Corn	16" x 24" Food Tray
No. 715B	7" x 15"	Soybean	17" x 19" Seedburo Germ.
No. 715C	7" x 15"	Corn	17" x 19" Seedburo Germ.



PURITY WORKBOARD AND DIAPHANOSCOPE

Seedburo's combination Purity Workboard and Diaphanoscope provides a convenient work area on a smooth opal glass surface for seed and grain analysis. A light placed under and shining through the glass assists in critical examination of the sample. Used by grain inspectors and seed analysts, the Purity Workboard and Diaphanoscope is ideal for grain grading and

checking for stress cracks (use fluorescent model for this application), germ damage, and seed viability. Unit features a removable sample drawer. Net wt. 12 lbs. Actual dims 32" L x 10½" W x 4¼" H. Ship wt 19 lbs. Ship dims 35" L x 17" W x 8" H.

Purity Workboard & Diaphanoscope with Incandescent Bulb, 110V, 50/60Hz	No. 135
Purity Workboard & Diaphanoscope with Incandescent Bulb, 220V, 50/60Hz	No. 135/A
Corn Stress Crack Viewer with Fluorescent Bulb, 110V, 60Hz	No. 135F/B
Replacement Opal Glass (only)	No. 135G



SEEDBURO HAND TALLY

This four-digit hand counter records to 9999. White numerals on a black background for easy reading. Thumb pressure on button produces audible click with each increase in count. Hand Tally features a large finger ring and turn-off (zero) knob. Allows for greater accuracy in seed counts and many other applications. Actual dims 1¾" dia. x 1-1/8"W.

Seedburo Hand Tally No. HT



DIGITAL HAND TALLY

A hand held counter with 4 digital LCD display and large positive action count button actuator. To reset all you need to do is hold the On/Off button for 3 seconds. A "beeper", which can be disabled, sounds every count to verify a count has been registered. Made of high impact resistant plastic combining light weight with outstanding durability. Operates on 1 AG10 1.5V battery (included) or equivalent. Dims. 2.4" L x 1.4" W x 0.5" D. Wt. 0.7 oz.

Digital Hand Tally No. HTD



SEEDBURO TETRAZOLIUM POWDER

Used by seed testing laboratories and universities; 2, 3 5 Triphenyl Tetrazolium Chloride is a seed testing compound used to determine seed viability. Distinguishing between viable and dead tissues of the embryo on the basis of their relative respiration rate in the hydrated state.

Tetrazolium Powder, 50 gram jar No. 3126



SEEDBURO WEED SEED BOOKLET

Ideal for lab or classroom, the Seedburo Weed Seed Booklet is completely illustrated with all of the common noxious and poisonous weeds and weed seeds with commonly used crop seeds.

Seedburo Weed Seed Booklet No. WSB